## NEDERLANDS

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Presentatie: 20 – 50 gram tubes

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**Presentation: 20 - 50 gram tubes**

### **Product description**

*Melloxy® contains:* honey, ozonated olive oil, propylene glycol and PEG 4000. When the gel is in contact with the wound, fluid is extracted from the surrounding tissues. The osmotic action creates a moist wound-healing environment. Honey together with the aforementioned moist wound healing environment stimulates the wound healing process, facilitates autolysis and promotes better epithelial cell migration. The gel will dilute gradually and can be removed easily.

Melloxy® is a primary wound dressing that can be covered by most commonly used secondary dressings.

### **Complete quantitative and qualitative composition**

Melloxy® consists of 40% honey. The composition below may differ slightly because honey is a natural product.

Nr.	component	%	min - max
1	Water	7,7	6,7 - 9
2	Fructose	17,5	12,5 - 18
3	Glucose	14	8,5 - 14,5
4	Maltose	3	2 - 4

5	Sucrose	1,8	0,05 - 2
6	PEG4000	16	0
7	Propylene glycol	34,5	0
8	Palmitic acid	1,08	0
9	Palmitoleic acid	0.06	0.05 – 0.065
10	Stearic acid	0.33	0.30 – 0.35
11	Oleic acid	5.3	5.4 – 5.5
12	Linoleic acid	0.73	0.63 – 0.8
13	Nonanal	0.36	0.3 – 0.4
14	Pelargonic acid	3.5	3.4 – 3.6
15	Azelaaldehydic Acid	0.36	0.31- 0.4
16	Capro aldehyde	0.19	0.14 – 0.23
17	Azelaic Acid	0.27	0.22 – 0.30
18	Non-3-enoic acid	0.3	0.25 – 0.35
19	Caproic acid	0.02	0.018 – 0.021
20	Undecylic acid	0.009	0.008 – 0.010

### Indications and Intended use

*Superficial Chronic wounds* - Decubitus, Ulcus Cruris, diabetic foot wounds, Fungating wounds (deodorising and debriding).

*Contaminated Acute wounds* – surgical wounds (postoperative wounds), Traumatic wounds (superficial wounds, cuts), Small burn and laser wounds (1° and 2° degree). The ointment creates a moist wound climate.

Melloxy® can be used both by professionals and laymen.

### Contraindications

No contraindications have been reported to date.

Do not use the gel if the patient is highly prone to infection e.g. a patient with very large burn wounds.

Do not use on patients with a known hypersensitivity to the gel or its components.

Do not use on eyes or mucous membranes.

Do not use if package is damaged.

### Side Effects

- Melloxy® can cause in some cases a short stinging sensation.

- If a serious incident has occurred with regard to the medical device (Melloxy®), the patient/user should report this to the manufacturer and the competent authority of the Member State where the user and/or patient.

### Precautions and observations

- Do not use if the tube has been damaged.

- Do not use the product continuously for more than 30 days.

- The presence of bacteria in wounds is a normal phenomenon and is no contra indication.

- If an infection is suspected or is developing during the use then a correct antibiotic therapy should be started. The use of Melloxy® can nevertheless be continued provided that the evolution of the wound is carefully monitored and that the treatment is placed under medical supervision.

- Depending on the underlying pathology the wound care treatment should always be accompanied by the correct medical treatment and supervision.

- Melloxy® creates a moist wound climate which stimulates the debridement of necrotic tissue. Due to this it is normal that the wound surface becomes bigger.

- When the condition of the wound deteriorates, treatment with Melloxy® should be stopped.

- We recommend that a tube is used on one patient to avoid cross-contamination (= single patient device).
- Always use Melloxy® exactly as your nurse, doctor or pharmacist has told you, or as explained in this insert.

### **Pregnancy, breastfeeding and use in Pediatric patients (0-15 years)**

There is no information available about the use of Melloxy® during pregnancy and breastfeeding. Are you pregnant or Breastfeeding? Please contact your doctor or pharmacist. No contraindications are known to date for the use of Melloxy® in pediatric patients.

### **Instructions for use**

1. Clean the wound.
2. Remove lid and open tube in an aseptic way.
3. Squeeze the tube to extract the gel and apply it without touching the wound.
4. Use a (sterile) spatula when necessary.
5. Place a dressing on top of the gel to cover the wound.
6. Melloxy® can remain on the wound up to 48 hours on average. Reapply the gel on the wound as soon as the gel has disappeared.

### **Remove**

Melloxy® should be changed when leakage of the secondary dressing occurs or when the secondary dressing needs to be removed. The ointment can be cleansed out of the wound with a sterile saline solution. Dressing change is depending on the wound fluid production of the wound, but is usual between 2 -7 days.

### **Waste processing/disposal**

Cardboard box, paper leaflet for paper waste and the plastic tube for recyclable waste.

### **Storage precautions**

- Store between 5°C - 25° C
- After opening the package, the product can be used for 28 days.

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© SanoMed Manufacturing BV 2022 Barring Spelling errors, SanoMed Manufacturing BV shall at all times retain the right to change its products and assortment.